




Instrument Name: ZenoTOF™ 7600 System LC200
Instrument Model: ZenoTOF 7600 System
Ion Source: OptiFlow 1-50µL Micro/MicroCal
Serial Number: FB24962403

Manufacturer: AB Sciex
Target Instrument: ZenoTOF 7600 System
Optimization Set: Negative Quick Status Check

1. Negative MS Check

Step	Result	Details																														
Achieve Stable Spray / Modify	Pass	Slope: 0.000490027127601987, previous was 0.000490027127601987 Delay: -13.5975402795579, previous was -13.5975402795579																														
Average of 13 spectra - Confirmation Mode Complete																																
																																
<table><tr><th>Item</th><th>Min</th><th>Max</th><th>Actual</th><th>Status</th></tr><tr><td>Mass (Da)</td><td>112.8856</td><td>113.0856</td><td>112.9845</td><td>Pass</td></tr><tr><td>Sum Intensity (cps)</td><td>1000</td><td>1e+11</td><td>18302</td><td>Pass</td></tr></table>	Item	Min	Max	Actual	Status	Mass (Da)	112.8856	113.0856	112.9845	Pass	Sum Intensity (cps)	1000	1e+11	18302	Pass	<table><tr><th>Item</th><th>Min</th><th>Max</th><th>Actual</th><th>Status</th></tr><tr><td>Mass (Da)</td><td>792.7596</td><td>792.9596</td><td>792.8526</td><td>Pass</td></tr><tr><td>Sum Intensity (cps)</td><td>1000</td><td>1e+11</td><td>10187</td><td>Pass</td></tr></table>	Item	Min	Max	Actual	Status	Mass (Da)	792.7596	792.9596	792.8526	Pass	Sum Intensity (cps)	1000	1e+11	10187	Pass	
Item	Min	Max	Actual	Status																												
Mass (Da)	112.8856	113.0856	112.9845	Pass																												
Sum Intensity (cps)	1000	1e+11	18302	Pass																												
Item	Min	Max	Actual	Status																												
Mass (Da)	792.7596	792.9596	792.8526	Pass																												
Sum Intensity (cps)	1000	1e+11	10187	Pass																												
																																

Channel Alignment

Pass

Sum of 10 spectra

Instrument Name: ZenoTOF™ 7600 System LC200
Instrument Model: ZenoTOF 7600 System
Ion Source: OptiFlow 1-50µL Micro/MicroCal
Serial Number: FB24962403


Manufacturer: AB Sciex
Target Instrument: ZenoTOF 7600 System
Optimization Set: Negative Quick Status Check

TOF MS Mass Check

Pass

Slope: 0.000490027727589735, previous was 0.000490027127601987
Delay: -13.6079799911381, previous was -13.5975402795579

Average of 10 spectra - Confirmation Mode Complete




TOF MS/MS Mass Check

Pass

Slope: 0.000490026677628662, previous was 0.000490025856063894
Delay: -13.7206011802538, previous was -13.6955712782684

Average of 10 spectra - Confirmation Mode Complete



TOF MS/MS Mass Check Zeno On

Pass

Slope: 0.00049002599877421, previous was 0.000490025589254436
Delay: -13.7051245473183, previous was -13.6839087004933

Average of 10 spectra - Confirmation Mode Complete

