

Ion Source:

## **Instrument Tuning Report**

**nstrument Name:** ZenoTOF™ 7600 System LC200

Instrument Model: ZenoTOF 7600 System

OptiFlow 1-50µL Micro/MicroCal

Serial Number: FB24962403

Manufacturer: AB Sciex

**Target Instrument:** ZenoTOF 7600 System

Optimization Set: Negative Quick Status Check

#### 1. Negative MS Check


#### Step Result Details

Achieve Stable Spray / Modify


Pass Slope: 0.000490027997223775, previous was

0.000490027997223775

Delay: -13.6609218388699, previous was -13.6609218388699



# Channel Alignment Pass Sum of 10 spectra 792.8596 792.8596



Align to channel:	z ke	Suit:	Pass		
Item	Min	Max	Initial	Actual	Status
Channel 1 vs. 2 Slope	-50	50	0.047056451	0.044144507	Pass
Channel 1 vs. 2 Intercept	-50	50	-10.756244091	-10.29417986	Pass
Channel 3 vs. 2 Slope	-50	50	-0.108683695	-0.115973246	Pass
Channel 3 vs. 2 Intercept	-50	50	-4.561174419	-3.700061918	Pass
Channel 4 vs. 2 Slope	-50	50	-0.286370769	-0.286891018	Pass
Channel 4 vs. 2 Intercept	-50	50	7.84049216	7.864162853	Pass

www.SCIEX.com Page 1 of 2



**Instrument Model:** 

TOF MS Mass Check

Ion Source:

### **Instrument Tuning Report**

ZenoTOF™ 7600 System LC200

ZenoTOF 7600 System

OptiFlow 1-50µL Micro/MicroCal

Serial Number: FB24962403

-2.00

Error (ppm)

2.00

-0.12554 Pass

-2.00

Error (ppm)

2.00

-1.33638 Pass

Manufacturer: AB Sciex


**Target Instrument:** ZenoTOF 7600 System

**Optimization Set:** Negative Quick Status Check

Slope: 0.000490027127601987, previous was **Pass** 


0.000490027997223775

Delay: -13.5975402795579, previous was -13.6609218388699



TOF MS/MS Mass Check **Pass** Slope: 0.000490025856063894, previous was 0.000490027252191172

Delay: -13.6955712782684, previous was -13.700717174369



TOF MS/MS Mass Check Zeno On **Pass** Slope: 0.000490025589254436, previous was 0.000490026696015616 Delay: -13.6839087004933, previous was -13.6742995868526 Average of 10 spectra - Confirmation Mode Complete 248.9604 384.9352 520.9100 656.8848 384.9347 520.9103-2500 200 300 2000 150 400 200 1500 100 1000 100 200 500 244 245 246 247 248 249 250 251 252 253 516 517 518 519 520 521 522 523 524 525 383 384 385 386 387 388 389 655 656 657 Item Min Max Actual Status Mass (Da) 248.8604 249.0604 248.9604 Pass Mass (Da) 384.8352 385.0352 384.9347 Pass Mass (Da) 520.8100 521.0100 520.9103 Pass Mass (Da) 656.7848 656.9848 656.8841 Pass 100000 30564 Pass 24000 Resolution 22000 Resolution 100000 30035 Pass Resolution 25000 100000 34208 Resolution 26000 100000 33931 Pass Pass Sum Intensity Sum Intensity Sum Intensity Sum Intensity 6233.2 Pass (cps) (cps) (cps)

Error (ppm)

-2.00

2.00

0.59593 Pass

-2.00

Error (ppm)

2.00

-1.02119 Pass

www.SCIEX.com Page 2 of 2