

Ion Source:

Instrument Tuning Report

Instrument Name: ZenoTOF™ 7600 System

Instrument Model: ZenoTOF 7600 System

OptiFlow 50-200µL Micro/MicroCal

Serial Number: FB24962403

Manufacturer: AB Sciex

Target Instrument: ZenoTOF 7600 System

Optimization Set: Negative Quick Status Check

1. Negative MS Check

Step Result Details


Achieve Stable Spray / Modify

Pass

Slope: 0.000489981816387641, previous was

0.000489981816387641

Delay: -13.5880122297166, previous was -13.5880122297166



Pass Channel Alignment Sum of 10 spectra 112.9856 792.8596 112.9846 792.8550 2e+05 3e+05 2e+05 1e+05 792.4 1126 112.8 113 113.2 113.4 792.6 m/z, Da m/z, Da Item Min Max Mass (Da) 112.7356 113.2356 112.9846 Mass (Da) 792.6096 793.1096 792.855 Resolution 0 100000 34909 Pass Resolution 0 100000 44385 Pass Sum Intensity Sum Intensity 1e+11 1.3277e+06 Pass 0 1e+11 1.9298e+06 Pass (cps) Align to channel: 2 Result: Pass Min Max Initial Actual Channel 1 vs. 2 Slope -50 50 0.037243019 0.039439446 Pass Channel 1 vs. 2 Intercept -50 50 -10.059662371 -10.325723086 Pass Channel 3 vs. 2 Slope -0.100170143 -0.098559883 Pass -4.237146866 -4.591876019 Pass Channel 4 vs. 2 Slope -0.260238406 -0.258449241 Channel 4 vs. 2 Intercept -50 50 7 572741476 7 095849877

www.SCIEX.com Page 1 of 2



Instrument Model:

TOF MS Mass Check

Ion Source:

Instrument Tuning Report

nstrument Name: ZenoTOF™ 7600 System

ZenoTOF 7600 System

OptiFlow 50-200µL Micro/MicroCal

Serial Number: FB24962403

Manufacturer: AB Sciex


Target Instrument: ZenoTOF 7600 System

Optimization Set: Negative Quick Status Check

Pass Slope: 0.000489982501683296, previous was

0.000489981816387641


Delay: -13.5888724244004, previous was -13.5880122297166



TOF MS/MS Mass Check Pass Slope: 0.000489981123283425, previous was

0.000489980320682435


Delay: -13.684879475171, previous was -13.6787985214466



TOF MS/MS Mass Check Zeno On Pass Slope: 0.000489980772376221, previous was

0.000489980033858822

Delay: -13.6576103675257, previous was -13.6540099212559



www.SCIEX.com Page 2 of 2