

Instrument Tuning Report

ZenoTOF™ 7600 System **Instrument Name:**

Instrument Model: ZenoTOF 7600 System

OptiFlow 50-200µL Micro/MicroCal Ion Source:

FB24962403 Serial Number:

Manufacturer: **AB Sciex**


Target Instrument: ZenoTOF 7600 System

Optimization Set: Positive Quick Status Check

1. Positive MS Check

Step Result **Details** Achieve Stable Spray / Modify

Pass Slope: 0.000490102425189461, previous was



Pass Channel Alignment Sum of 10 spectra 132.9049 829.5393 2e+05 132.9044-829.5378 2e+05 ფ 1.5e+05 සු 1.5e+05 ntensity, 1e+05 5e+04 5e+04 132.8 133 133.2 133.4 829 2 829.4 829.6 Item Min Max Actual Mass (Da) 132.6549 133.1549 132.9044 Mass (Da) 829.2893 829.7893 829.5378 Resolution 0 100000 36852 Pass Resolution 0 100000 40325 Pass Sum Intensity Sum Intensity 1e+11 7.851e+05 Pass 0 1e+11 1.9637e+06 Pass (cps) Align to channel: 2 Result: Pass Min Max Initial Actual Channel 1 vs. 2 Slope -50 50 -0.126227253 -0.125456582 Pass Channel 1 vs. 2 Intercept -50 50 -10.542220504 -10.753156757 Pass Channel 3 vs. 2 Slope -3.631603363 -3.749719316 Pass Channel 4 vs. 2 Slope 0.116171189 0.1195739 Channel 4 vs. 2 Intercept -50 50 9 272105403 8 702500278

www.SCIEX.com Page 1 of 2



Instrument Model:

Instrument Tuning Report

ZenoTOF™ 7600 System

ZenoTOF 7600 System

Target Instrument: ZenoTOF 7600 System

OptiFlow 50-200µL Micro/MicroCal Ion Source: **Optimization Set:** Positive Quick Status Check

Serial Number: FB24962403


TOF MS Mass Check **Pass** Slope: 0.00049010248950908, previous was

AB Sciex

0.000490102425189461

Manufacturer:


Delay: -13.4517952035879, previous was -13.46034855723



TOF MS/MS Mass Check **Pass** Slope: 0.000490101551950563, previous was

0.000490101678289791


Delay: -13.4808324288948, previous was -13.4777224599844



TOF MS/MS Mass Check Zeno On Slope: 0.000490102775985005, previous was **Pass**

0.000490102702566802

Delay: -13.464094087671, previous was -13.4605135351005



www.SCIEX.com Page 2 of 2