

Ion Source:

## **Instrument Tuning Report**

**nstrument Name:** ZenoTOF™ 7600 System Optiflow

Instrument Model: ZenoTOF 7600 System

OptiFlow 50-200µL Micro/MicroCal

Serial Number: FB24962403

Manufacturer: AB Sciex

Target Instrument:ZenoTOF 7600 System

**Optimization Set:** Positive Quick Status Check

## 1. Positive MS Check


 Step
 Result
 Details

 Achieve Stable Spray / Modify
 Pass
 Slope: 0

Pass Slope: 0.000490048306655696, previous was

0.000490048306655696

Delay: -13.3814662735105, previous was -13.3814662735105



## **Pass Channel Alignment** Sum of 10 spectra 132.9049 829.5393 132.9014 2.5e+05 829.522 නු 1.5e+05 2e+05 1.5e+05 Intensity, 1e+05 1e+05 5e+04 132.8 133 133.2 133.4 829 2 829.4 829.6 Item Max Actual Mass (Da) 132.6549 133.1549 132.9014 Mass (Da) 829.2893 829.7893 829.5227 Resolution 0 100000 39567 Pass Resolution 0 100000 42049 Pass Sum Intensity Sum Intensity 1e+11 7.289e+05 Pass 1e+11 2.1852e+06 Pass (cps) Align to channel: 2 Result: Pass Min Max Initial Actual Status Channel 1 vs. 2 Slope -50 50 -0.121138847 -0.114577356 Pass Channel 1 vs. 2 Intercept -50 50 -11.282765973 -11.520622 Pass Channel 3 vs. 2 Slope 0.071488354 0.075846596 Pass -3.882108725 -4.660597023 Pass Channel 4 vs. 2 Slope 0.108195137 0.111743367 Channel 4 vs. 2 Intercept -50 50 8.403234367 8 686778703

www.SCIEX.com Page 1 of 2



**Instrument Model:** 

TOF MS Mass Check

## **Instrument Tuning Report**

**Pass** 

ZenoTOF™ 7600 System Optiflow

ZenoTOF 7600 System

OptiFlow 50-200µL Micro/MicroCal Ion Source:

Serial Number: FB24962403 Manufacturer: AB Sciex


ZenoTOF 7600 System **Target Instrument:** 

Positive Quick Status Check **Optimization Set:** 

Slope: 0.000490052117855052, previous was

0.000490048306655696


Delay: -13.4509500231973, previous was -13.3814662735105



TOF MS/MS Mass Check **Pass** Slope: 0.000490051338487275, previous was

0.000490047561614517


Delay: -13.4729870054648, previous was -13.4458611459108



TOF MS/MS Mass Check Zeno On **Pass** Slope: 0.000490052717189903, previous was

0.00049004865812057

Delay: -13.4217190111928, previous was -13.4355774428237



www.SCIEX.com Page 2 of 2