

Ion Source:

Instrument Tuning Report

ZenoTOF™ 7600 System

Instrument Model: ZenoTOF 7600 System

OptiFlow 50-200µL Micro/MicroCal

FB24962403 Serial Number:

Manufacturer: **AB Sciex**

Target Instrument: ZenoTOF 7600 System

Optimization Set: Positive Quick Status Check

1. Positive MS Check


Step Result **Details Pass**

Achieve Stable Spray / Modify

Slope: 0.00049011915417191, previous was

0.00049011915417191

Delay: -13.4211588177695, previous was -13.4211588177695



Pass Channel Alignment Sum of 10 spectra 132.9049 829.5393 132,9027 1.5e+05 1.5e+05 1e+05 1e+05 5e+04 5e+04 132.8 133 133.2 133.4 829 2 829.4 829.6 Item Min Max Actual Mass (Da) 132.6549 133.1549 132.9027 Mass (Da) 829.2893 829.7893 829.5297 Resolution 0 100000 33962 Pass Resolution 0 100000 36478 Pass Sum Intensity Sum Intensity 1e+11 7.8993e+05 Pass 0 1e+11 1.9333e+06 Pass (cps) Align to channel: 2 Result: Pass Min Max Initial Actual Channel 1 vs. 2 Slope -50 50 -0.126298304 -0.125020072 Pass Channel 1 vs. 2 Intercept -50 50 -11.413154947 -10.208557454 Pass Channel 3 vs. 2 Slope 0.076260159 0.076531736 Pass -3.752435641 -3.755834018 Pass Channel 4 vs. 2 Slope 0.131201677 0.12869308 Channel 4 vs. 2 Intercept -50 50 8 429155041 9.084837768

www.SCIEX.com Page 1 of 2



Instrument Model:

TOF MS Mass Check

Ion Source:

Instrument Tuning Report

Pass

ZenoTOF™ 7600 System

ZenoTOF 7600 System

OptiFlow 50-200µL Micro/MicroCal **Optimization Set:** Positive Quick Status Check

Manufacturer:

Serial Number: FB24962403


ZenoTOF 7600 System **Target Instrument:**

AB Sciex

Slope: 0.000490121080626083, previous was

0.00049011915417191

Delay: -13.4666123420202, previous was -13.4211588177695



TOF MS/MS Mass Check **Pass** Slope: 0.000490120020005401, previous was 0.000490118201860901 Delay: -13.493605766298, previous was -13.4980239786222 Average of 10 spectra - Confirmation Mode Complete 185.1285 494.3337 607.4178 1.2e+04 185.1285 298.2127 494.3341 607.4179 1200 5000 1e+04 8000 1000 4000 8000 6000 800 3000 600 611.4127 2000 4000 400 498.3286 1000 2000 200 187.1077 183 184 185 186 188 294 295 296 297 298 490 491 492 493 494 603 604 605 606 607 608 609 610 611 299 300 301 Min Max Actual Item Min Max Actual Status Item Max Actual Status Item Status Item Min Max Actual Status Pass Mass (Da) 185 0285 185 2285 185 1285 Pass Mass (Da) 298 1125 298 3125 298 2127 Pass Mass (Da) 494.2337 494.4337 494.3341 Mass (Da) 607.3178 607.5178 607.4179 Pass Resolution 21000 100000 23306 Pass Resolution 22000 100000 29384 Pass Resolution 25000 100000 31697 Pass Resolution 26000 100000 30457 Pass Sum Intensity Sum Intensity Sum Intensity Sum Intensity 100 1e+11 9690.6 Pass 100 1e+11 40615 Pass 100 1e+11 99077 Pass 100 1e+11 92392 Pass Error (ppm) -2 00 2.00 0.21505 Pass Error (ppm) -2 00 2.00 0.73998 Pass Error (nnm) -2 00 2.00 0.74249 Pass Error (ppm) -2 00 2.00 0.15921 Pass

TOF MS/MS Mass Check Zeno On **Pass** Slope: 0.000490121030609039, previous was 0.000490119449824009 Delay: -13.4776990437519, previous was -13.4903094474059 Average of 10 spectra - Confirmation Mode Complete 185.1285 298.2125 494.3337 607.4178 494.3337--2500 3000 2000 600 2000 1500 2000 1500 400 611.4121 1000 1000 1000 200 500 498 3283 500 0 183 184 186 187 188 189 190 294 295 296 297 298 490 491 492 493 494 495 496 497 603 604 605 606 607 608 609 610 611 612 185 299 300 Item Max Actual Item Max Item Min Max Max Mass (Da) 185.0285 185.2285 185.1283 Pass Mass (Da) 298.1125 298.3125 298.2124 Pass Mass (Da) 494.2337 494.4337 494.3337 Pass Mass (Da) 607.3178 607.5178 607.4174 Pass Resolution 21000 100000 27137 Pass Resolution 22000 100000 31093 Pass Resolution 25000 100000 30313 Pass Resolution 26000 100000 31379 Sum Intensity Sum Intensity Sum Intensity Sum Intensity 28377 100 23433 (cps) (cps) (cps) (cps) Error (ppm) -2.00 2.00 -0.88991 Pass Error (ppm) -2.00 2.00 -0.44303 Pass Error (ppm) -2.00 2.00 -0.09143 Pass Error (ppm) -2.00 2.00 -0.61316 Pass

www.SCIEX.com Page 2 of 2