

**Instrument Model:** 

Ion Source:

## **Instrument Tuning Report**

**Pass** 

**Instrument Name:** ZenoTOF™ 7600 System LC200

ZenoTOF 7600 System

OptiFlow 1-50µL Micro/MicroCal

Serial Number: FB24962403

Manufacturer: AB Sciex

**Target Instrument:** ZenoTOF 7600 System

Optimization Set: Negative Quick Status Check

### 1. Negative MS Check


Step Result Details

Achieve Stable Spray / Modify

Slope: 0.000489990657945629, previous was


0.000489990657945629

Delay: -13.5955573632082, previous was -13.5955573632082



### Channel Alignment

#### **Pass**



Item	Min	Max	Initial	Actual	Status
Channel 1 vs. 2 Slope	-50	50	0.036181878	0.041129906	Pass
Channel 1 vs. 2 Intercept	-50	50	-10.061376086	-10.43249889	Pass
Channel 3 vs. 2 Slope	-50	50	-0.107695509	-0.107366206	Pass
Channel 3 vs. 2 Intercept	-50	50	-3.585505072	-3.94961802	Pass
Channel 4 vs. 2 Slope	-50	50	-0.274855849	-0.269955971	Pass
Channel 4 vs. 2 Intercept	-50	50	8.355386412	7.910273767	Pass

www.SCIEX.com Page 1 of 2



**Instrument Model:** 

TOF MS Mass Check

# **Instrument Tuning Report**

**nstrument Name:** ZenoTOF™ 7600 System LC200

ZenoTOF 7600 System

**Ion Source:** OptiFlow 1-50µL Micro/MicroCal

Serial Number: FB24962403

Manufacturer: AB Sciex


Target Instrument: ZenoTOF 7600 System

**Optimization Set:** Negative Quick Status Check

Pass Slope: 0.000489990275006645, previous was

0.000489990657945629


Delay: -13.5789137258397, previous was -13.5955573632082




TOF MS/MS Mass Check Pass Slope: 0.000489989015638218, previous was

0.000489989453612945

Delay: -13.6832899700256, previous was -13.6794256221511



Delay: -13.6751130110889, previous was -13.6410189235612



www.SCIEX.com Page 2 of 2