



**Instrument Name:** ZenoTOF™ 7600 System  
**Instrument Model:** ZenoTOF 7600 System  
**Ion Source:** OptiFlow 50-200µL Micro/MicroCal  
**Serial Number:** FB24962403


**Manufacturer:** AB Sciex  
**Target Instrument:** ZenoTOF 7600 System  
**Optimization Set:** Negative Quick Status Check

## TOF MS Mass Check

**Pass**

Slope: 0.000490006061144034, previous was 0.000490007055250039  
Delay: -13.5917812909534, previous was -13.5923578824717

### Average of 10 spectra - Confirmation Mode Complete




## TOF MS/MS Mass Check

**Pass**

Slope: 0.000490004976083629, previous was 0.000490005797420302  
Delay: -13.6708964026432, previous was -13.6807293681015

### Average of 10 spectra - Confirmation Mode Complete



## TOF MS/MS Mass Check Zeno On

**Pass**

Slope: 0.000490004611772901, previous was 0.000490005377863692  
Delay: -13.634454984552, previous was -13.6760491653842

### Average of 10 spectra - Confirmation Mode Complete

