



**Instrument Name:** ZenoTOF™ 7600 System  
**Instrument Model:** ZenoTOF 7600 System  
**Ion Source:** OptiFlow 50-200µL Micro/MicroCal  
**Serial Number:** FB24962403


**Manufacturer:** AB Sciex  
**Target Instrument:** ZenoTOF 7600 System  
**Optimization Set:** Positive Quick Status Check

## TOF MS Mass Check

**Pass**

Slope: 0.000490065877906141, previous was 0.000490063407953435  
Delay: -13.4378937856746, previous was -13.3984681442952

### Average of 10 spectra - Confirmation Mode Complete




## TOF MS/MS Mass Check

**Pass**

Slope: 0.000490064757364887, previous was 0.000490062548042322  
Delay: -13.4780633380219, previous was -13.4722631610112

### Average of 10 spectra - Confirmation Mode Complete



## TOF MS/MS Mass Check Zeno On

**Pass**

Slope: 0.000490065734569162, previous was 0.000490063988451274  
Delay: -13.4504216509354, previous was -13.4436180227862

### Average of 10 spectra - Confirmation Mode Complete

