




**Instrument Name:** ZenoTOF™ 7600 System  
**Instrument Model:** ZenoTOF 7600 System  
**Ion Source:** OptiFlow 50-200µL Micro/MicroCal  
**Serial Number:** FB24962403

**Manufacturer:** AB Sciex  
**Target Instrument:** ZenoTOF 7600 System  
**Optimization Set:** Negative Quick Status Check

## 1. Negative MS Check

Step	Result	Details																														
Achieve Stable Spray / Modify	Pass	Slope: 0.000490109621396227, previous was 0.000490109621396227 Delay: -13.6443087348115, previous was -13.6443087348115																														
Average of 30 spectra - Confirmation Mode Complete																																
																																
<table><tr><th>Item</th><th>Min</th><th>Max</th><th>Actual</th><th>Status</th></tr><tr><td>Mass (Da)</td><td>112.8856</td><td>113.0856</td><td>112.9843</td><td>Pass</td></tr><tr><td>Sum Intensity (cps)</td><td>1000</td><td>1e+11</td><td>7362.1</td><td>Pass</td></tr></table>	Item	Min	Max	Actual	Status	Mass (Da)	112.8856	113.0856	112.9843	Pass	Sum Intensity (cps)	1000	1e+11	7362.1	Pass	<table><tr><th>Item</th><th>Min</th><th>Max</th><th>Actual</th><th>Status</th></tr><tr><td>Mass (Da)</td><td>792.7596</td><td>792.9596</td><td>792.8526</td><td>Pass</td></tr><tr><td>Sum Intensity (cps)</td><td>1000</td><td>1e+11</td><td>3119.9</td><td>Pass</td></tr></table>		Item	Min	Max	Actual	Status	Mass (Da)	792.7596	792.9596	792.8526	Pass	Sum Intensity (cps)	1000	1e+11	3119.9	Pass
Item	Min	Max	Actual	Status																												
Mass (Da)	112.8856	113.0856	112.9843	Pass																												
Sum Intensity (cps)	1000	1e+11	7362.1	Pass																												
Item	Min	Max	Actual	Status																												
Mass (Da)	792.7596	792.9596	792.8526	Pass																												
Sum Intensity (cps)	1000	1e+11	3119.9	Pass																												
																																

## Channel Alignment

Pass

Sum of 10 spectra

</

**Instrument Name:** ZenoTOF™ 7600 System  
**Instrument Model:** ZenoTOF 7600 System  
**Ion Source:** OptiFlow 50-200µL Micro/MicroCal  
**Serial Number:** FB24962403


**Manufacturer:** AB Sciex  
**Target Instrument:** ZenoTOF 7600 System  
**Optimization Set:** Negative Quick Status Check

## TOF MS Mass Check

**Pass**

Slope: 0.000490111265480814, previous was 0.000490109621396227  
Delay: -13.6724345418674, previous was -13.6443087348115

### Average of 10 spectra - Confirmation Mode Complete




## TOF MS/MS Mass Check

**Pass**

Slope: 0.000490110439986943, previous was 0.000490108628261658  
Delay: -13.6926500657664, previous was -13.7104336672888

### Average of 10 spectra - Confirmation Mode Complete



## TOF MS/MS Mass Check Zeno On

**Pass**

Slope: 0.000490110228058364, previous was 0.000490108427821228  
Delay: -13.6652994621755, previous was -13.6748317007687

### Average of 10 spectra - Confirmation Mode Complete

