




Instrument Name: ZenoTOF™ 7600 System
Instrument Model: ZenoTOF 7600 System
Ion Source: OptiFlow 50-200µL Micro/MicroCal
Serial Number: FB24962403

Manufacturer: AB Sciex
Target Instrument: ZenoTOF 7600 System
Optimization Set: Positive Quick Status Check

1. Positive MS Check


Step	Result	Details																													
Achieve Stable Spray / Modify	Pass	Slope: 0.000490184163622215, previous was 0.000490184163622215 Delay: -13.4677977697352, previous was -13.4677977697352																													
Average of 30 spectra - Confirmation Mode Complete																															
																															
<table><tr><th>Item</th><th>Min</th><th>Max</th><th>Actual</th><th>Status</th></tr><tr><td>Mass (Da)</td><td>132.8049</td><td>133.0049</td><td>132.9057</td><td>Pass</td></tr><tr><td>Sum Intensity (cps)</td><td>1000</td><td>1e+11</td><td>10402</td><td>Pass</td></tr></table>	Item	Min	Max	Actual	Status	Mass (Da)	132.8049	133.0049	132.9057	Pass	Sum Intensity (cps)	1000	1e+11	10402	Pass	<table><tr><th>Item</th><th>Min</th><th>Max</th><th>Actual</th><th>Status</th></tr><tr><td>Mass (Da)</td><td>829.4393</td><td>829.6393</td><td>829.5442</td><td>Pass</td></tr><tr><td>Sum Intensity (cps)</td><td>1000</td><td>1e+11</td><td>17007</td><td>Pass</td></tr></table>	Item	Min	Max	Actual	Status	Mass (Da)	829.4393	829.6393	829.5442	Pass	Sum Intensity (cps)	1000	1e+11	17007	Pass
Item	Min	Max	Actual	Status																											
Mass (Da)	132.8049	133.0049	132.9057	Pass																											
Sum Intensity (cps)	1000	1e+11	10402	Pass																											
Item	Min	Max	Actual	Status																											
Mass (Da)	829.4393	829.6393	829.5442	Pass																											
Sum Intensity (cps)	1000	1e+11	17007	Pass																											
																															

Channel Alignment

Pass


Sum of 10 spectra

132.9049



Item	Min	Max	Actual	Status
Mass (Da)	132.6549	133.1549	132.9052	Pass
Resolution	0	100000	38056	Pass
Sum Intensity (cps)	0	1e+11	6.6443e+05	Pass

829.5393



Item	Min	Max	Actual	Status
Mass (Da)	829.2893	829.7893	829.5431	Pass
Resolution	0	100000	42993	Pass
Sum Intensity (cps)	0	1e+11	1.8779e+06	Pass

Align to channel: 2 Result: Pass

Item	Min	Max	Initial	Actual	Status
Channel 1 vs. 2 Slope	-50	50	-0.119882282	-0.119170367	Pass
Channel 1 vs. 2 Intercept	-50	50	-11.649330296	-11.296168377	Pass
Channel 3 vs. 2 Slope	-50	50	0.073929345	0.076649551	Pass
Channel 3 vs. 2 Intercept	-50	50	-3.66568469	-4.065933879	Pass
Channel 4 vs. 2 Slope	-50	50	0.119244883	0.119755202	Pass
Channel 4 vs. 2 Intercept	-50	50	8.705899509	8.926913551	Pass

Instrument Name: ZenoTOF™ 7600 System
Instrument Model: ZenoTOF 7600 System
Ion Source: OptiFlow 50-200µL Micro/MicroCal
Serial Number: FB24962403


Manufacturer: AB Sciex
Target Instrument: ZenoTOF 7600 System
Optimization Set: Positive Quick Status Check

TOF MS Mass Check

Pass

Slope: 0.000490182509969467, previous was 0.000490184163622215
Delay: -13.4649883170905, previous was -13.4677977697352

Average of 10 spectra - Confirmation Mode Complete




TOF MS/MS Mass Check

Pass

Slope: 0.000490181366635546, previous was 0.000490183224885963
Delay: -13.4980766626877, previous was -13.4953577985652

Average of 10 spectra - Confirmation Mode Complete



TOF MS/MS Mass Check Zeno On

Pass

Slope: 0.000490182874592014, previous was 0.000490184964817994
Delay: -13.4565460981884, previous was -13.4706142660205

Average of 10 spectra - Confirmation Mode Complete

